II<sup>a</sup> Classe su filo ortodontico  
II<sup>a</sup> Class on orthodontic wire

II<sup>a</sup> Classe su bande  
II<sup>a</sup> Class on bands

III<sup>a</sup> Classe su bande  
III<sup>a</sup> Class on bands

## Il Second e il Third Class Resolver (S&TCR)

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*S&TCR are two fixed orthopaedic devices functioning with proprioceptive stimulus, using the same principles as other orthopaedic devices, but thanks to their 24 hours-a-day action, they get the therapeutical result independently from the patient cooperation, in shorter time. The two devices act as functional systems because of their work, simultaneous to the stomatognathic apparatus functions (mastication, deglutition, phonation).*

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
### *S&TCR for the latero-deviation solution*

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*S&TCR, in particular the SCR version, can be easily used with a successful result in the latero-deviation Second Class resolution, acting through an asymmetrical activation. Working 24 hours-a-day, this device corrects the latero-deviation problem thanks to the continuous stimulus coming from the muscular receptors that have the maximum possibility to change the neuromuscular engram, helped by the neuronal tissue plasticity.*

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### *Use of SCR on .021 x .025 to correct a lateral deviation*

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*SCR version is highly indicated for the resolution of second class latero-deviation, through an asymmetrical activation and working with multibrackets devices when setted on .021 x .025 or .019 x .025 wires. Working 24 hours-a-day, this device corrects the latero-deviation taking advantage of the continuous stimulus from the muscular receptors which has, thanks to the neuronal tissue plasticity, the maximum possibility to modify the engram.*

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


## Second Class Resolver (SCR)

The Second Class Resolver (SCR) does not force the mandible to keep the same position (as the Herbst does), because the arches stay independent of each other and very important factor: the stimulus inducted to the mandible is proprioceptive. Thanks to the proprioceptive contact between the two sliding planes of the device (upper Skibar and lower Blockrun) and to the fact that the condyle and the masseters are displaced of about 2 mm, the SCR succeeds in orientating either the muscular or the condylar growth vectors. The therapeutical result concretizes the mandibular branch advancement and the growth of the mandibular ante-rotation.

- The device results very comfortable for the patient; after 4-7 days the patient adapts to it easily.
- The device in its 3 different variants (on multibrackets, sectional 4-6 on the 4 quadrants, buccal splint from 6 to 6, upper and lower) is placed in the patient vestibule, leaving the oral cavity free.
- The oral cavity free is advantageous for mastication, deglutition and phonation, even when the logopaedic rehabilitative therapy or the physiotherapy are requested. This situation represents the innovation of this device: the possibility of using an orthopaedic device for the correction of the skeletal class and at the same time the performance of myofunctional and/or logopaedic therapies. All perfectly compatible with the use of orthodontic appliances (brackets, bands, etc. ). The activation of SCR is easy: it is enough to act only on the tightening screw to mobilize the Blockrun. The Blockrun distal sliding activates the device. The brand new way of interfacing of both planes determines the mandibular advancement.

## Third Class Resolver (TCR)

The Third Class Resolver is a fixed orthopaedic device which does not force the mandible to keep the same position, to be used in growth age. It uses the same orthopaedic devices principles as the stimulation of the circum-maxillary sutures with the aim of promoting the upper mandibular advancement. The characteristics of this device are:

1. Working 24 hours-a-day it enables the therapeutic result independently from the patient cooperation, in shorter time.
2. Thanks to the TCR design, the stimuli, promoting the mandibular advancement, are activated by the act of deglutition; for this reason, probably, they are characterized by light forces;
3. The device biomechanics enables the dynamic which transmits the force to the maxilla to stimulate its advancement. This action/reaction force is mostly confined inside the device. This is summarized by a slight reaction retrusive force, loading on the mandible. Thanks to the proprioceptive contact between the two sliding planes of the device, the TCR enables the orientation of the muscular and the condylar growth vectors. Considering that the main body of the device S&TCR does not invade the oral cavity, the TCR is one of the most valid therapeutical instruments to resolve the skeletal Class III in patients in growth age, even contemporarily at a myofunctional therapy. This situation portrays the innovation of the device: the possibility of using an orthopaedic device for the correction of the skeletal class and at the same time the performance of myofunctional and/or logopaedic therapies. All perfectly compatible with the use of orthodontic appliances (brackets, bands, etc. ). Stimulating, then, for the orthodontist who can chose several therapeutical combinations, in order to satisfy all patients needs, shortening treatment times. The device results comfortable for the patient; after 4-7 days the patient adapts to it easily. The device, in its 2 variants (sectional 4-6 on the 4 quadrants; buccal splint from 6 to 6, upper and lower) is placed behind the coilspring, which is placed behind the Blockrun, in the lower part of the device. The mesial sliding of the stop enables the whole system advancement activating the device.