

EASY BURN




DENTAL WORLD
ITALIAN MANUFACTURER

Brucia aghi e strumenti acuminati
Sharpened instruments and needle burner device

Il Brucia Aghi Easy Burn, in grado di distruggere aghi e piccoli strumenti acuminati e taglienti, per evitare il rischio di infezione. Un piccolo cratere con una temperatura di 1400°C, per incenerire all'istante l'ago appena usato e con esso ogni possibilità di contaminazione. La batteria consente un'autonomia di 400 - 600 aghi prima della ricarica; con caricabatterie esterno che garantisce ogni sicurezza di normativa CE. I residui combusti da fusione, trovano posto nel contenitore interno, facilmente ispezionabile per la pulizia.

Easyburn is a needle burner capable of small needles and sharp instruments destruction, avoiding infection risks. A small crater with a 1400°C temperature provides the necessary heat capable of instant needle destruction, avoiding any contamination. The battery has a duration of 400-600 needles burns before the recharge process; the device has an external CE regulation - compliant battery charger. The garbages of the burning process go to an internal easily attainable cleanable container.



Brucia Aghi Easy Burn 00340710



Dental World s.r.l. Via delle Ricamatrici B/4 Zona PIP, 70056 Molfetta (BA) ITALY

tel. +39 080.338.10.04 / fax +39 080.338.66.17 - www.dental-world.it / e.mail: info@dental-world.it