



Bande Molari GX / Molar Bands

Le bande anatomiche Ortho World per molari e premolari, sono realizzate con il migliore acciaio medicale. Disponibili nelle misure da 1 a 32 per i molari mascellari e mandibolari e da 1 a 24 per i premolari (universali). La marchiatura a laser permette una facile identificazione delle misure mentre, l'interno sabbiato, permette un'ottima ritenzione ai cementi ortodontici. Garantiamo inoltre un'eccellente servizio di prepuntatura.

Anatomical bicuspid and molar bands, made with the best medical steel. Size 1 to 32 for molar bands; 1 to 24 for bicuspid bands (universal). Permanent laser identification markings and sandblasted for accurate cement adhesion. Dental World offers precision prewelding of bands.

Assortement Kit

band size	01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Trial kit 80 pcs	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	4	2	-	-	-	-	-	-	-	-	-	-	-	-
Intro kit 100 pcs	2	2	2	2	2	4	4	4	4	4	6	6	6	4	4	4	4	4	4	2	2	2	2	2	2	2	2	2	2	2	2	
Pro Kit.150 pcs	2	2	2	2	4	4	6	6	8	8	8	10	10	10	8	8	8	6	6	4	4	2	2	2	2	2	2	2	2	2	2	
Master Kit 300 pcs	4	4	6	6	8	8	14	14	16	16	18	18	18	18	16	16	14	14	8	8	6	6	6	4	4	4	4	4	4	4	4	

I Kit sono formati da 50% di bande destre e 50% di bande sinistre / 50 % right bands and 50% left bands

Molar Bands Kit

	Trial Kit 80 pcs.	Intro Kit 100 pcs.	Pro Kit 150 pcs.	Master Kit 300 pcs.
UPPER	KMBU180	KMBU1100	KMBU1150	KMBU1300
LOWER	KMBL180	KMBL1100	KMBL1150	KMBL1300

Bande Anatomiche GX / Anatomical bands



UL	MBUL1*
UR	MBUR1*

* aggiungere sempre la misura desiderata Es. MBUL1 25

* insert desired molar band size Ex. MBUL1 25



LL	MBLL1*
LR	MBLR1*

* aggiungere sempre la misura desiderata Es. MBLL1 25

* insert desired molar band size Ex. MBLL1 25

Bande Prepuntate GX / Prewelded bands



ROTH	Torque -10°	Offset 8°	.022"
UL	MBUL1* + WBT2210LG		
UR	MBUR1* + WBT2210RG		

* aggiungere sempre la misura desiderata Es. MBUL1 25

* insert desired molar band size Ex. MBUL1 25



ROTH	Torque -25°	Offset 8°	.022"
LL	MBLL1* + WBT2225LG		
LR	MBLR1* + WBT2225RG		

* aggiungere sempre la misura desiderata Es. MBLL1 25

* insert desired molar band size Ex. MBLL1 25



EDGEWISE	Torque 0°	Offset 0°	.022"
UL	MBUL1* + WBT2200LG		
UR	MBUR1* + WBT2200RG		

* aggiungere sempre la misura desiderata Es. MBUL1 25

* insert desired molar band size Ex. MBUL1 25



EDGEWISE	Torque 0°	Offset 0°	.022"
LL	MBLL1* + WBT2200LG		
LR	MBLR1* + WBT2200RG		

* aggiungere sempre la misura desiderata Es. MBLL1 25

* insert desired molar band size Ex. MBLL1 25



RICKETTS	Torque -10°	Offset 8°	.022"
UL	MBUL1* + WBT3810LG		
UR	MBUR1* + WBT3810RG		

* aggiungere sempre la misura desiderata Es. MBUL1 25

* insert desired molar band size Ex. MBUL1 25



RICKETTS	Torque 0°	Offset 0°	.022"
LL	MBLL1* + WBT2800L		
LR	MBLR1* + WBT2800R		

* aggiungere sempre la misura desiderata Es. MBLL1 25

* insert desired molar band size Ex. MBLL1 25

■ Merce disponibile / Goods available

■ Merce a richiesta / Goods on request

Bande Premolari GX / Bicuspid Bands

Assortement Kit

UPPER		01	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
band size																									
Intro kit 100 pcs		2	2	2	2	2	4	4	4	6	6	8	8	8	8	6	6	4	4	4	2	2	2	2	2
Pro Kit.150 pcs		2	2	2	4	4	6	6	6	8	8	10	10	10	10	10	8	8	8	6	4	4	2	2	

LOWER		03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
band size																									
Intro kit 100 pcs		2	2	2	2	2	4	4	4	6	6	8	8	8	8	6	6	4	4	4	2	2	2	2	2
Pro Kit.150 pcs		2	2	2	4	4	6	6	6	8	8	10	10	10	10	10	8	8	8	6	4	4	2	2	

I Kit sono formati da 50% di bande destre e 50% di bande sinistre / 50% right bands and 50% left bands

Bicuspid Bands Kit

	Intro Kit 100 pcs.	Pro Kit 150 pcs.
UPPER (size 1-24)	OUBBK100	OUBBK150
LOWER (size 3-26)	OLBBK100	OLBBK150

Ricambio 5pz.
OUBB*
OLBB*

* aggiungere sempre la misura desiderata Es. OUBB 14

* insert desired premolar band size Ex. OUBB 14



Ortho World Tubi da puntare non convertibili / Non-convertible weldable tubes

I Tubi sono prodotti con tecnologia MIM di altissima qualità in confezioni da 10 pezzi. Possono essere forniti singolarmente o saldati sulle bande.

High quality and precision tubes realized with MIM technology. They can be purchased loose or pre-welded to our bands.



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	WBT1800R	WBT1800L	WBT1200R	WBT1200L
ROTH	-10°	8°	WBT1810R	WBT1810L	WBT1210R	WBT1210L
ROTH	-25°	8°	WBT1825R	WBT1825L	WBT1225R	WBT1225L
10 pcs / pack						



Double .045 Occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	WBT2800RO	WBT2800LO	WBT2200RO	WBT2200LO
ROTH	-10°	8°	WBT2810RO	WBT2810LO	WBT2210RO	WBT2210LO
10 pcs / pack						



Double .045 Gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	WBT2800RG	WBT2800LG	WBT2200RG	WBT2200LG
ROTH	-10°	8°	WBT2810RG	WBT2810LG	WBT2210RG	WBT2210LG
ROTH	-25°	8°	WBT2825LG	WBT2825RG	WBT2225RG	WBT2225LG
10 pcs / pack						



Double Bioform			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
RICKETTS	0°	0°	WBT2800R	WBT2800L	-	-
10 pcs / pack						



Triple .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
RICKETTS	-10°	8°	WBT3810RG	WBT3810LG	-	-
10 pcs / pack						

Ortho World Tubi da puntare convertibili / *Convertible weldable tubes*



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
ROTH	-10°	8°	WBT1810RC	WBT1810LC	WBT1210RC	WBT1210LC
ROTH	-25°	8°	WBT1825LC	WBT1825RC	WBT1225LC	WBT1225RC
10 pcs / pack						



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
MBT	-20°	0°	CS1-816	CS1-815	CS1-216	CS1-215
10 pcs / pack						



Double auxiliary			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
ROTH	-25°	8°	WBT2825LC	WBT2825RC	WBT2225LC	WBT2225RC
10 pcs / pack						



Double .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	8°	WBT2800RCO	WBT2800LCO	WBT2200RCO	WBT2200LCO
ROTH	-10°	5°	WBT2810RCO	WBT2810LCO	WBT2210RCO	WBT2210LCO
10 pcs / pack						



Double .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
ROTH	-25°	8°	WBT2825LCG	WBT2825RCG	WBT2225LCG	WBT2225RCG
10 pcs / pack						



Doppio .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
MBT	-14°	10°	CC1-817	CC1-818	CC1-217	CC1-218
10 pcs / pack						



Triple .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
MBT	-14°	10°	CT1-807	CT1-808	CT1-207	CT1-208
10 pcs / pack						



Triple .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
ROTH	-10°	8°	WBT3810RCO	WBT3810LCO	WBT3210RCO	WBT3210LCO
10 pcs / pack						

Ortho World Tubi da cementare non convertibili / *Bondable non-convertible tubes*

Prodotto di altissima qualità realizzato con tecnologia MIM, bassetta ritentiva e con disegno anatomico, per adattarsi in maniera ottimale ai molari.

High quality product, manufactured with MIM technology. The bonding base are anatomically designed and contoured for positive adaptation to molars.



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	8°	BBT18008R	BBT18008L	BBT12008R	BBT12008L
ROTH	-10°	8°	BBT18108R	BBT18108L	BBT12108R	BBT12108L
ROTH	-25°	5°	BBT18255R	BBT18255L	BBT12255R	BBT12255L
10 pcs / pack						



Double .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	8°	BBT28008OR	BBT28008OL	BBT22008OR	BBT22008OL
ROTH	-10°	5°	BBT28108OR	BBT28108OL	BBT22108OR	BBT22108OL
10 pcs / pack						



Double Bioform			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
RICKETTS	0°	0°	BBT2800R	BBT2800L	-	-
10 pcs / pack						

Ortho World Tubi da cementare convertibili / *Bondable convertible tubes*



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
ROTH	-10°	8°	BBT18108RC	BBT18108LC	BBT12108RC	BBT12108LC
ROTH	-25°	5°	BBT18255LC	BBT18255RC	BBT12255LC	BBT12255RC
10 pcs / pack						



Double			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
ROTH	-25°	5°	BBT28255LC	BBT28255RC	BBT22255LC	BBT22255RC
10 pcs / pack						



Double .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
ROTH	-10°	8°	BBT28108ORC	BBT28108OLC	BBT22108ORC	BBT22108OLC
10 pcs / pack						



Double .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
ROTH	-10°	8°	BBT28108GRC	BBT28108GLC	BBT22108GRC	BBT22108GLC
10 pcs / pack						



Double .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
ROTH	-25°	5°	BBT28255GLC	BBT28255GRC	BBT22255GLC	BBT22255GRC
10 pcs / pack						

Tubi mini da cementare 1° & 2° molare / 1° & 2° molar bondable mini tube



Single non convertible			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	61-711-11	61-711-21	61-712-11	61-712-21
ROTH - MAXILLARY	-10°	8°	61-721-11	61-721-21	61-722-11	61-722-21
ROTH - MANDIBULAR	-25°	5°	61-721-31	61-721-41	61-722-31	61-722-41
MBT - MAXILLARY	-14°	10°	61-731-11	61-731-21	61-732-11	61-732-21
MBT - MANDIBULAR	-20°	0°	61-731-31	61-731-41	61-732-31	61-732-41
10 pcs / pack						

Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	NS2B-821	NS2B-822	NS2B-221	NS2B-222
10 pcs / pack						



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
ROTH	-10°	0°	NS2B-831	NS2B-832	NS2B-231	NS2B-232
10 pcs / pack						

Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
MBT	-14°	0°	NS2B-835	NS2B-836	NS2B-235	NS2B-236
10 pcs / pack						

Single			.018x.025		.022x.028	
	TORQUE	OFFSET	LL	LR	LL	LR
ROTH / MBT	-25°	0°	NS2B-834	NS2B-833	NS2B-234	NS2B-233
10 pcs / pack						

BioLine Tubi da cementare nickel-free / Nickel free bondable tubes

I tubi BioLine Nickel-Free rappresentano l'unica alternativa per pazienti sensibili al Nickel.

Nickel-Free bondable tubes represent the alternative product for patients who experience sensitivity to Nickel products.



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
MBT	-14°	8°	BBTM1814R	BBTM1814L	BBTM1214R	BBTM1214L
MBT	-20°	0°	BBTM1820L	BBTM1820R	BBTM1220L	BBTM1220R
5 pcs / pack						



Single			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
ROTH	-14°	10°	BBT1401	BBT1402	BBT1451	BBT1452
ROTH	-25°	4°	BBT1414	BBT1413	BBT1464	BBT1463
5 pcs / pack						

SENSITIVE
NICKEL-FREE
PATIENTS

■ Merce disponibile / Goods available

■ Merce a richiesta / Goods on request

Bande Molari TX / Molar Bands

Le bande anatomiche TX per molari, sono realizzate con il migliore acciaio medicale. Disponibili nelle misure da 29 a 44 per i molari mascellari e mandibolari. La marchiatura a laser permette una facile identificazione delle misure mentre, l'interno sabbiato, permette un'ottima ritenzione ai cementi ortodontici. Garantiamo inoltre un'eccellente servizio di prepuntatura.

Molar anatomical bands, made with the best medical steel. Sizes 29 to 44. Permanent laser identification markings and sandblasted for accurate cement adhesion. Dental World offers precision prewelding of bands.

Assortement Kit

band size	29.5	30	30.5	31	31.5	32	32.5	33	33.5	34	34.5	35	35.5	36	36.5	37	37.5	38	38.5	39	39.5	40	40.5	41	41.5	42	42.5	43	43.5	44
Trial kit 80 pcs	-	-	-	-	-	2	2	4	4	4	6	6	10	10	10	6	6	4	4	2	-	-	-	-	-	-	-	-	-	-
Pro Kit.150 pcs	2	2	2	2	2	2	4	4	6	8	10	10	10	12	12	12	10	8	6	4	4	4	2	2	2	2	2	2	2	2
Master Kit 300 pcs	8	8	8	8	8	8	16	16	24	32	40	40	40	40	48	48	40	32	24	16	16	16	8	8	8	8	8	8	8	8

Molar Bands Kit

	Trial Kit 80 pcs.	Pro Kit 150 pcs.	Master Kit 300 pcs.
UPPER	KMBEU180	KMBEU1150	KMBEU1300
LOWER	KMBEL180	KMBEL1150	KMBEL1300

Bande Anatomiche TX / Anatomical bands



UL	MBEUL1*
UR	MBEUR1*

* aggiungere sempre la misura desiderata Es. MBEUL1 29.5
* insert desired molar band size Ex. MBEUL1 29.5



LL	MBELL1*
LR	MBELR1*

* aggiungere sempre la misura desiderata Es. MBLL1 29.5
* insert desired molar band size Ex. MBLL1 29.5

Bande Prepuntate TX / Prewelded bands

Le bande TX possono essere prepuntate a scelta con i tubi sotto elencati, nelle varie tecniche. Per i kit e i ricambi, basta aggiungere alla fine del codice il riferimento del tubo scelto.

(Es. MBEUL129 + 64-611-12)

banda + doppio tubo

Bands can be prewelded with different technique tubes. For kit and refills, add tube reference at the end of code

(Es. MBEUL129 + 64-611-12)

band + double tube



ROTH	Torque -10°	Offset 8°	.022"
UL	MBEUL1* + 64-622-22		
UR	MBEUR1* + 64-622-12		

* aggiungere sempre la misura desiderata Es. MBEUL1 29.5
* insert desired molar band size Ex. MBEUL1 29.5



ROTH	Torque -25°	Offset 5°	.022"
LL	MBELL1* + 64-622-32		
LR	MBELR1* + 64-622-42		

* aggiungere sempre la misura desiderata Es. MBELL1 29.5
* insert desired molar band size Ex. MBELL1 29.5



MBT	Torque -14°	Offset 10°	.022"
UL	MBEUL1* + 64-622-22		
UR	MBEUR1* + 64-622-12		

* aggiungere sempre la misura desiderata Es. MBEUL1 29.5
* insert desired molar band size Ex. MBEUL1 29.5



MBT	Torque -20°	Offset 0°	.022"
LL	MBELL1* + 64-632-32		
LR	MBELR1* + 64-632-42		

* aggiungere sempre la misura desiderata Es. MBELL1 29.5
* insert desired molar band size Ex. MBELL1 29.5

- Merce disponibile / Goods available
- Merce a richiesta / Goods on request

Eclipse Tubi da puntare / *Weldable tubes*



Single non convertible			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	63-611-11	63-611-21	63-612-11	63-612-21
ROTH - MAXILLARY	-10°	8°	63-621-11	63-621-21	63-622-11	63-622-21
ROTH - MANDIBULAR	-25°	5°	63-621-31	63-621-41	63-622-31	63-622-41
MBT - MAXILLARY	-14°	10°	63-631-11	63-631-21	63-632-11	63-632-21
MBT - MANDIBULAR	-20°	0°	63-631-31	63-631-41	63-632-31	63-632-41
10 pcs / pack						



Single convertible			.018x.025		.022x.028	
	TORQUE	OFFSET	UR	UL	UR	UL
ROTH - MAXILLARY	-10°	8°	64-621-11	64-621-21	64-622-11	64-622-21
ROTH - MANDIBULAR	-25°	5°	64-621-41	64-621-31	64-622-41	64-622-31
MBT - MAXILLARY	-14°	10°	64-631-11	64-631-21	64-632-11	64-632-21
MBT - MANDIBULAR	-20°	0°	64-631-41	64-631-31	64-632-41	64-632-31
10 pcs / pack						



Double convertible .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	64-611-12	64-611-22	64-612-12	64-612-22
ROTH - MAXILLARY	-10°	8°	64-621-12	64-621-22	64-622-12	64-622-22
ROTH - MANDIBULAR	-25°	5°	64-621-32	64-621-42	64-622-32	64-622-42
MBT - MAXILLARY	-14°	10°	64-631-12	64-631-22	64-632-12	64-632-22
MBT - MANDIBULAR	-20°	0°	64-631-32	64-631-42	64-632-32	64-632-42
10 pcs / pack						



Triple convertible .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
ROTH - MAXILLARY	-10°	8°	64-621-13	64-621-23	64-622-13	64-622-23
MBT - MAXILLARY	-14°	0°	64-631-13	64-631-23	64-632-13	64-632-23
10 pcs / pack						

Eclipse Tubi da cementare / *Bondable tubes*



Single non convertible			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	61-611-11	61-611-21	61-612-11	61-612-21
ROTH - MAXILLARY	-10°	8°	61-621-11	61-621-21	61-622-11	61-622-21
ROTH - MANDIBULAR	-25°	5°	61-621-31	61-621-41	61-622-31	61-622-41
MBT - MAXILLARY	-14°	10°	61-631-11	61-631-21	61-632-11	61-632-21
MBT - MANDIBULAR	-20°	0°	61-631-31	61-631-41	61-632-31	61-632-41
10 pcs / pack						



Double convertible .045 gingival			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
EDGEWISE	0°	0°	62-611-12	62-611-22	62-612-12	62-612-22
ROTH - MAXILLARY	-10°	8°	62-621-12	62-621-22	62-622-12	62-622-22
ROTH - MANDIBULAR	-25°	5°	62-621-32	62-621-42	62-622-32	62-622-42
MBT - MAXILLARY	-14°	10°	62-631-12	62-631-22	62-632-12	62-632-22
MBT - MANDIBULAR	-20°	0°	62-631-32	62-631-42	62-632-32	62-632-42
10 pcs / pack						



Triple convertible .045 occlusal			.018x.025		.022x.028	
	TORQUE	OFFSET	UR/LL	UL/LR	UR/LL	UL/LR
ROTH - MAXILLARY	-10°	8°	62-621-13	62-621-23	62-622-13	62-622-23
MBT - MAXILLARY	-14°	0°	62-631-13	62-631-23	62-632-13	62-632-23
10 pcs / pack						

Accessori per bande / Band attachments

Lingual Sheaths

Realizzati per accettare Barre Palatali, Quad Helix o altri format removibili con diametri fino a .036" x .072". Devono essere assicurati con una legatura metallica o elastica per evitare che vengano ingeriti.

Designed to accept palatal bars, Quad Helix appliances and other removable appliances incorporating .036" x .072" distal engagement ends. To be locked with metal or elastomeric ligature ties.



WLS055R

WLS055L





WLS155

W/Hook

Lingual Sheaths left	WLS055L
Lingual Sheaths right	WLS055R

W/OHook

Universal	WLS155
10 pcs / pack	



Lingual Cleat

Ideali per agganciare legature metalliche o elastiche.

Ideal for metal and elastomeric ligatures.



Cleat linguale short	WLCSHORT
Cleat linguale long	WLCLONG
10 pcs / pack	

Brackets da puntare / Weldable brackets

Brackets da puntare su bande evadibili nelle prescrizioni ROTH ed EDGEWISE.

Weldable Brackets - Roth and Edgwise technique.

Bicuspid Brackets

ROTH

	TORQUE	OFFSET	.018"	.022"
Upper	-7°	0°	WRRXU518	WRRXU522
Lower	-17°	0°	WRRXL418	WRRXL422
Lower, 2nd Bicuspid	-22°	10°	WRRXL518	WRRXL522
10 pcs / pack				

STANDARD

	TORQUE	OFFSET	.018"	.022"
Universal	0°	0°	WEDBR18	WEDBR22
10 pcs / pack				

Molar Brackets

STANDARD EDGEWISE

	TORQUE	OFFSET	.018"	.022"
Brackets Convertible	0°	0°	EDSWMBC018	EDSWMBC022
Brackets Non-Convertible	0°	0°	EDSWMB018	EDSWMB022
10 pcs / pack				



PUNTATRICE ORTODONTICA


- PUNTATURA / WELDING

- BRASATURA / BRAZE WELDING

- PUNTATURA A FREQUENZA / FREQUENCY WELDING

- RISCALDAMENTO E TEMPRA DEI FILI / WIRE HEAT TREATMENT

Bande e Tubi / Bands and Tube



Merce disponibile / Goods available

Merce a richiesta / Goods on request