



acusjet

Ago dentale per anestesia - Dental needle for anesthesia
MADE IN ITALY


MADE IN ITALY





acusjet

Ago dentale per anestesia - Dental needle for anesthesia

Descrizione / Characteristics

Acus Jet, aghi sottoposti a controlli qualitativi molto severi in ogni singola fase di produzione; presentano una cannula in acciaio inossidabile a parete ultrasottile con punta a triplice affilatura siliconata. Acus Jet sono aghi con cono filettato in polipropilene, sterilizzati a raggi gamma, con lotto di produzione, scadenza e misura stampati direttamente sull'ago per una veloce identificazione.

Acus Jet, needles exposed to very strict quality controls in each step of their production; they are made up of an ultra thin and stainless steel canula with a triple silicone sharpening tip. Acus Jet needles have a polypropylene threaded cone, gamma ray sterilized, with production line batch, expiry date and size printed directly on the needle for a faster identification.

Vantaggi / Advantages

L'indicatore di orientamento della punta permette al professionista di conoscere l'esatta direzione di uscita dell'anestetico.

The orientation indicator of the tip allows the professional to know the right exit direction of the anesthetic.

Packaging
















- Scatola: 100 aghi

Box: 100 needles

- Cartone: 12 scatole

Carton: 12 boxes

Codici aghi / Needle codes

0,4 x 8 mm		Extra corto - intraligamentare <i>Extra short - intraligamentary</i>		00200006	Acus Jet 27G - 0,4x8mm - 100pcs - Violet
0,3 x 12 mm		Corto - intraligamentare <i>Short - intraligamentary</i>		00200002	Acus Jet 30G - 0,3x12mm - 100pcs - White
0,3 x 16 mm		Corto - intraligamentare <i>Short - intraligamentary</i>		00200003	Acus Jet 30G - 0,3x16mm - 100pcs - Bordeaux
0,4 x 16 mm				00200008	Acus Jet 27G - 0,4x16mm - 100pcs - Pink
0,3 x 21 mm		Lungo - Infiltrazione <i>Long - infiltration</i>		00200004	Acus Jet 30G - 0,3x21mm - 100pcs - Blue
0,4 x 21 mm				00200009	Acus Jet 27G - 0,4x21mm - 100pcs - Yellow
0,3 x 25 mm		Lungo - Infiltrazione <i>Long - infiltration</i>		00200005	Acus Jet 30G - 0,3x25mm - 100pcs - Green
0,4 x 25 mm				00200010	Acus Jet 27G - 0,4x25mm - 100pcs - Orange
0,4 x 38 mm		Lungo - mandibolare <i>Long - mandibular</i>		00200012	Acus Jet 27G - 0,4x38mm - 100pcs - Light Blue