


PRO ENDO Micro Mirror

The Charm of *reflection*

Ami
ENDOSURGERY



The individual mirrors are interchangeable




Pro Endo Micro Mirror Dental World: innovative instruments designed to ensure precision, practicality, and the highest quality during dental procedures.

KEY FEATURES

- Micro mirrors with optimized sizes ranging from 2 to 9 mm, available in rectangular and round shapes.
- Ergonomic design with a non-slip handle.
- Durable and easy to replace interchangeable tips.
- High-quality materials that ensure durability and reliability.
- Special coating: cold-plated.
- Anti-drip and anti-fog surface.
- Scratch-resistant.
- Autoclavable at 134°C.



Perfect reflection
even when in contact



Shapeable tip

KIT - 4 pcs

790080-KIT Pro Endo Micro Mirror Kit

MIRRORS - 1 pc

790080-C02	Mirror C02 - Round Ø3mm
790080-C04	Mirror C04 - Round Ø5mm
790080-R02	Mirror R02 - Square 2x7mm
790080-R04	Mirror R04 - Square 3x9mm

MIRRORS WITH HANDLES - 1 pc

790080-C02M	Mirror with handles C02M - Round Ø3mm
790080-C04M	Mirror with handles C04M - Round Ø5mm
790080-R02M	Mirror with handles R02M - Square 2x7mm
790080-R04M	Mirror with handles R04M - Square 3x9mm