

FILLETCH

Mordenzante smalto dentinale 37% - **SENZA SILICE**
 37% enamel-dentin etching agent - **SILICA FREE**



BEFORE

AFTER

precise etching.
 controlled application.

- Semi gel (a base di polimeri naturali)
Gel seeds (made from natural polymers)
- Tissotropico
Thixotropic
- Nessuna ostruzione dell'agetto
No clogging of the tip
- Biocompatibile
Biocompatible
- Privo di silice
Silica-free



Descrizione - Description

Il mordenzante smalto dentinale al 37% senza silice Fill Etch è stato sviluppato per garantire un'adesione ottimale nei trattamenti clinici quotidiani. La sua speciale formulazione priva di silice rappresenta un vantaggio cruciale: garantisce una maggiore forza adesiva, favorendo una profonda ed efficace penetrazione degli adesivi direttamente nei tubuli dentinali. È specificamente indicato per la mordenzatura di smalto e dentina, così come per la mordenzatura di cementi vetroionomerici. Fill Etch si presenta sotto forma di un semi gel a base di polimeri naturali, dotato di un'eccellente natura tissotropica che permette un'applicazione precisa e controllata esattamente dove serve, senza il rischio di colature. Essendo un materiale biocompatibile e idrosolubile al 100%, assicura una completa detersione del campo operatorio non lasciando alcun residuo dopo la fase di lavaggio.

Fill Etch has been developed to ensure optimal adhesion in daily clinical treatments. Its special silica-free formulation represents a crucial advantage: it guarantees greater bond strength by promoting deep and effective penetration of adhesives directly into the dentinal tubules. It is specifically indicated for etching enamel and dentine, as well as for etching glass ionomer cements. Fill Etch comes as a semi-gel based on natural polymers, featuring an excellent thixotropic nature that allows for precise and controlled application exactly where needed, without the risk of running. As a 100% biocompatible and water-soluble material, it ensures complete cleansing of the operative field, leaving no residue after the rinsing phase.



Vantaggi - Benefits

- Idrosolubile al 100%
100% water-soluble
- Facile da applicare e rimuovere
Easy to apply and remove

APPLICATION



30ml x 1 syringe FILL Etch



5ml x 1 syringe FILL Etch

3ml x 2 syringes FILL Etch

3ml x 5 syringes FILL Etch



0000166	FILL Etch - 1 syringe x 5ml
0000168	FILL Etch - 2 syringes x 3ml
0000169	FILL Etch - 5 syringes x 3ml
0000170	FILL Etch - 1 syringe x 30ml